

DNase Footprint
(From Kathy Jones)

- 1) DNA Mix
 - 32P labeled probe
 - 10 ul 10% polyvinyl alcohol
 - 1 ul 10 A 260 unit/ml poly d(I-C) (from Sigma P9514)
 - H₂O to 25 ul
- 2) Add Z' or TM buffer to a final volume of 50 ul allowing room for the amount of protein extract used.
- 3) Add protein extract last.
- 4) incubate on ice 15 min.
- 5) DNase digestion. (amount of DNase to use must be titrated for best results)
 - a) 1 ul of 2.5 mg/ml DNase I is added to 1 ml ice cold water, mix thoroughly by inversion.
 - b) Add 100 ul from above to 200 ul ice cold water.
 - c) Add 50 ul 10 mM Mg Cl₂, 5mM Ca Cl₂ to the DNA protein mix leave at room temp for 1 min.
 - d) Add 2.0 ul diluted DNase, mix quickly , incubate at room temp. for 1 min. (If done at 1 min intervals, one can do three reactions at 15-20 second intervals.)
 - e) Add 90 ul STOP Mix, vortex, and store on ice.
 - f) After all reactions are completed to all;
 - Phenol extract.
 - Chloroform extract. This may require 1 min microfugation for good seperation.
 - Ethanol precipitate.
 - 70% ethanol wash and vacuum dry.
 - g) Add 4-10 ul loading buffer, boil for three minutes and load on sequencing gel.

Footprinting solutions.

	Stock	1 ml
TM buffer		
50mM Tris-HCl pH 7.9	2M	25 ul
12.5mM Mg Cl ₂	1M	12.5 ul
1mM EDTA	0.5M	2 ul

1mM DTT	1M	1 ul
20% glycerol	100%	200 ul
0.1M K Cl	1M	100 ul
H ₂ O	660 ul	

STOP Mix		5 ml
200mM Na Cl	1M	1 ml
20mM EDTA	0.5 M	250 ul
1% SDS	10%	500 ul
250 ug/ml tRNA	5 mg/ml	250 ul