

Tips for Growing Single Stranded DNA for Sequencing.

- 1) Maintain bacterial strains (JM101, JM103...) on minimal agarose plates.
- 2) Grow bacteria for as brief a time as possible.
- 3) Grow phage preps for 6-7 hours.
- 4) Do not freeze phage supernatents.
- 5) Unlike lambda, do not add chloroform to phage to store.

Minimal Media & Agarose

for 1 liter

- 15 gm. agarose [omit for media]
- 10.5 gm. K₂PO₄
- 4.5 gm. KH₂PO₄
- 1.0 gm [NH₄]SO₄
- 0.5 gm Na Citrate
- 0.2 gm MgSO₄ [sterilize separately/add when solution has cooled
some/ add 1 ml 1M MgSO₄]
- 5 mg thiamine HCl [filter sterilize separately]
- 2.0 mg glucose [sterilize separately/add 20 ml 20% glucose]

M13 Single Strand Template Purification for Sequencing.

1. Pick plaque and grow in 2.5ml 2xYT with 1-2 drops of fresh JM101 [JM103 or JM109 are also acceptable] for 6-7 hrs.
2. Spin out cells in microfuge 5 min. 2x.
3. To 1ml supernatant add 250ul 2.5M NaCl 20% PEG. Incubate 30 min. at R.T.
4. Spin in microfuge 5 min-pour off supernatant.
(Better- suck off with drawn out pasture pipette)
5. Quickly spin tubes again and remove residual supernatant.
6. Resuspend pellet in 100-200ml TES.
7. Extract: a) 2+ times phenol/chloroform mixture
b) chloroform
8. Add: a) 50ul Na Acetate
b) 250-500ul Etoh

- c) ppt. on dry ice.
- d) gently wash pellet with 70% Etoh.

9. Spin 5min. pour off supernatent & dry in vacuum.

10. Resuspend in 15 -20ul water.

Check 5ul on minigel.

Buffers for M13 purification.

20% Polyethanylene Glycol

2.5m NaCl [20gm PEG/100ml 2.5M NaCl]

TES

500mM NaCl	5M NaCl	50ml
10mM Tris pH7	2M Tris	2.5ml
1mM EDTA	.5M EDTA	1ml
	H2O	424ml
=====		
	Total	500ml

These solutions are used for preparing M13 single strand DNA for sequencing.