

## Dideoxy Sequencing Reactions

This is the original M13 sequencing protocol. All reagents must be prepared and calibrated before use. Recipes are provided at the end.

### 1) Annealing Sequencing Primer

- A) 5 ul DNA single stranded approx. 1ug
  - B) 1 ul Primer (17 bp)(NEB and BRL make suitable primers)
  - C) 1 ul 10xHin buffer
  - D) 5.5 ul H<sub>2</sub>O
- 12.5 ul Total Reaction

Heat to at least 65° for 5 min. (Not greater than 80° )  
Slow cool to room temp.

### 2) Add to each tube:

- 1 ul DTT 0.1 M
- 1 ul <sup>35</sup>S dATP
- 1 ul Klenow Pol I ( 1 unit per ul )
- 15.5 ul Total

### 3) Dilutions of deoxy nuclotides.

First Dilutions:

STOCK	DILUTION
dGTP (10mM)	4 --> 80
dATP (10mM)	4 --> 80
dTTP (10mM)	4 --> 80
dCTP (10mM)	4 --> 80
ddGTP (10mM)	1 --> 40
ddATP (10mM)	1 --> 160
ddTTP (10mM)	1 --> 10
ddCTP (10mM)	1 -->80

Second Dilutions ( mixes ):

	<u>dGTP</u>	<u>dATP</u>	<u>dTTP</u>	<u>dCTP</u>	<u>10 x HIN</u>
G°	1-2	-	20	20	20

A <sup>o</sup>	20	-	20	20	20
T <sup>o</sup>	20	-	1	20	20
C <sup>o</sup>	20	-	20	1	20
Chase	20	20	20	20	20

#### 4) Sequencing Reactions:

Set up 90 well conical bottom microtiter plate as follows;

(samples)	1	2	3...n		add 1ul	add 1 ul diluted
				G	G <sup>o</sup>	ddGTP
				A	A <sup>o</sup>	ddATP
			T	T <sup>o</sup>		ddTTP
		C	C <sup>o</sup>			ddCTP

add 3 ul of mix from step 2 to proper sample well.

add nucleotides to rows with corresponding names.

Incubate at 32<sup>o</sup> for 30 minutes. (This reaction can be run in wide range of temperatures from room temp. to 37<sup>o</sup>. There is also a protocol for sequencing DNA with G-C rich regions that uses a temp. of 42<sup>o</sup>.)

#### 5) Chase reaction:

a) add 1 ul chase (GATC)

b) 32<sup>o</sup> for 15 minutes.

#### 6) Stop reaction.

After chase add 5 ul gel dye and boil for 1 minute.

Ready to load on gel or freeze for later use.

#### DIDEOXY SEQUENCING BUFFERS

##### TM BUFFER

100mM TRIS HCl pH8.1	1M Tris HCl Ph8.1	100ul
100mM MgCl <sub>2</sub>	1M MgCl <sub>2</sub>	100ul
	H <sub>2</sub> O	800ul

	Total	1000ul
--	-------	--------

10x HIN BUFFER

70mM TRIS HCl pH7.6	2M Tris pH7.6	35ul
70mM MgCl <sub>2</sub>	1M MgCl <sub>2</sub>	70ul
500mM NaCl	5M NaCl	100ul
	H <sub>2</sub> O	795ul
	Total	1000ul

Note: Either buffer may be used in the priming reaction of M13 sequencing but TM seems to work slightly better.

Dideoxy or Maxam & Gilbert Sequencing Sample Dye Stop Buffer

0.5ml Formamide [deionized]

20ul 10% Bromphenol Blue

20ul 10% Xylene Cyanol

10ul 0.5ml EDTA

10xSTBE

	1 liter	2 Liter	4 Liter	6 Liter
Tris	108gm	216g	432gm	648gm
Boric Acid	55gm	110gm	220gm	330gm
EDTA	9.3gm	18.6gm	37.2g	55.8mg

Note: 10xSTBE-Sanger TBE used in gradient sequencing gels. This recipe is about twice the concentration of normal TBE.