

## Isolation of Ovarian Epithelial Cells

1. Remove ovaries from females
2. Label \_\_\_ eppendorfs
3. Add 300  $\mu$ l DMEM (w/o serum) to each
4. Add 3  $\mu$ l DiI-18\* to all except the controls for each
5. Place an ovary in each tube
6. Short time point: Tubes \_\_\_\_\_
  - leave for \_\_\_ minutes in the 37°C incubator
  - fix the ovary from tube \_\_\_ and the control with 4%PFA for \_\_\_ hrs. (Eppendorf)
  - transfer ovary from tube \_\_\_ to 100  $\mu$ l DMEM
    - triturate
    - fill the rest of the tube with 4%PFA
7. Long time point: Tubes \_\_\_\_\_
  - leave for \_\_\_ hours in the 37°C incubator
  - fix the ovary from tube \_\_\_ and the control with 4%PFA for \_\_\_ hrs. (Eppendorf)
  - transfer ovary from tube \_\_\_ to 100  $\mu$ l DMEM
    - triturate
    - fill the rest of the tube with 4%PFA
8. After each time point:
  - Rinse triturated pieces and whole ovaries well with 1X PBS
  - Cut the whole ovaries in half on a flat surface
  - Return them to the PBS
  - Take one half to histology
  - Use the other for vibratome
9. Use a transfer pipette to transfer the triturated pieces to a slide
  - Remove excess PBS
  - If pieces are too big, cut with a razor blade
  - Add a couple of drops of TOPO3: Vectashield mounting media (1:500)  
(Can also use 9 parts glycerol to 1 part PBS as mounting media)
  - View under confocal

\*dioctadecyl-6,6'-di(4-sulfophenyl)-3,3,3',3'-tetramethylindocarbocyanine  
(Molecular Probes)

## Vibratome

1. 7% - 9% agarose (similar consistency to tissue)
  - Microwave – turn off as soon as it starts to boil
  - Set in ~38°C - 43°C water bath for 15 min.
2. Vibratome supplies:
  - Ice bucket
  - Chilled 1X PBS
  - Blade (cut in half)
  - Microslide (shelf above vibratome)
  - Cover slip forceps (flat head)
  - Vice adaptor
  - Baster
  - Kimwipes
  - Superglue (Can get at storeroom – IRU 200)
3. Lower the stage and tighten using knob facing you
4. Clean a Microslide
5. Label the molding block w/ sharpie to keep track of your ovaries
  - Keep the same orientation throughout
6. Transfer tissues from PBS to the microslide
7. Fill molding block w/ agarose ~3/4 full
8. Working quickly
  - place tissue in agarose
  - orient so the blade hits the narrower side first
  - leave on ice for ~1 minute to solidify
9. Slide out the solidified agarose
10. Trim excess agarose w/ blade
11. Fill the black tray w/ ice
12. Make sure blade angle is ~23-25°
  - avg. speed: 2-3
  - avg. amp.: 8-9
13. Superglue the small piece of agarose containing your tissue onto the metal tray
14. Place in front of fan to dry (only for ~5 min.)
15. Place tray in slot and tighten using knob on the right
16. Fill with chilled PBS
17. Advance (FFW) until you reach the tissue
  - slow down when you get close to the tissue
  - adjust speed/amp as needed
18. Keep the specimens separate in the tray as they're slivered off
19. Store in PBS or transfer directly to a microscope slide
  - use TOPO3:Mounting Media (1:500)

\*dioctadecyl-6,6'-di(4-sulfohenyl)-3,3,3',3'-tetramethylindocarbocyanine  
(Molecular Probes)