

Isolation of Calvarial Osteoblasts

Procedure:

- 1) Cut out calvaria from P4 or P5 pups and place into CM (DMEM + 10 % FBS + 1 % Pen/Strep (P/S))
 - a. Euthanize pup by isoflourane inhalation
 - b. Douse pup with 70% EtOH
 - c. Collect tails for genotyping
 - d. Decapitate pup and cut the skin away
 - e. Hold the head with curved forceps placed through the orbita, cut the calvaria loose along the edge and place in a dish with PBS (+ 1 % P/S)
 - i. Cut away edges and sutures with scissors
 - f. Transfer calvaria to well of a 6-well dish containing PBS (+ 1 % P/S)
 - g. Repeat steps a – f for all remaining pups
 - h. Replace PBS with CM and incubate calvaria 1-2 days in 37 °C/ 5 % CO₂
- 2) Digest calvaria:
 - a. Make digest mix:
 - i. per genotype (~5 calvaria):
 1. 3 mg collagenase II
 2. 1 mL trypsin (0.25 % trypsin-EDTA; Gibco)
 3. up to 5.0 ml w/ PBS + 1 % P/S)...make 3X (per digest)
 4. **filter sterilize w/ 0.22 um filter...keep on ice**

number of genotype total =	3			
	Per genotype (ul)	Final conc per genotype	x (3.1) per digest	x total number of genotype (ml)
Collagenase (100mg/ml stock)	60	6 mg	186	0.558
Trypsin (.25%)	1000	~0.05%	3100	9.3
PBS + 1 % P/S	3940		12214	36.6

5 ml total

- b. Based on genotype, pool calvaria into 50 mL conical tube containing 5 mL PBS + P/S
 - i. Rinse 3X
- c. Add 4-5 mL digest mix to calvaria
 - i. Tighten cap and incubate at 37 °C for 5 min
 - ii. Discard supernatant
- d. Add fresh 5 mL digest mix to calvaria
 - i. With loose cap, incubate at 37 °C for 15 min
 - ii. Shake by hand for a few seconds
 - iii. Incubate at 37 °C for an additional 15 min
 - iv. Inactivate with 5 mL CM**
 - v. Save supernatant into new 50 ml conical tube (Incubate at 37 °C until step f is done)
- e. Add fresh 5 mL digest mix to calvaria
 - i. With loose cap, incubate at 37 °C for 15 min
 - ii. Shake by hand for a few seconds
 - iii. Incubate at 37 °C for an additional 15 min
 - iv. Inactivate with 5 mL CM**
 - v. Combine supernatant with that from step d.v.
- f. Spin down at 3000 rpm/5 min
- g. Aspirate supernatant
- h. Resuspend in DMEM + 10 % FBS + 1 % P/S or osteoblast differentiation media
 - i. Plate onto a well of a 6-well plate or a 60 mm dish (depending on size of cell pellet)
 - ii. Change media the following day
 - iii. Expand to 100% confluency for 1-2 weeks
 1. Feed every 3-4 days